

# HYDROGAMMA HCT

*energy-giving  
revitalizing  
detoxifying*

An ill-balanced diet, often combined with stress and lack of physical exercise, is one of the most frequent causes of problems related to the digestive system.

Some foods to blame for poor digestion/assimilation turn into "intestinal toxic waste" and this waste, over the course of time, ferments and release toxins which are then re-absorbed by the body, causing some of the most frequent diseases afflicting human kind (infections, polyarthrosis, rheumatism, psoriasis, hypertension...)

HYDROGAMMA HCT is a colon hydrotherapy system by means of which filtered water is infused into the colon via the rectum at body temperature, thereby obtaining a thorough, beneficial and complete cleansing action, which is both painless and odourless. The approach to this therapy is extremely gentle, both physically and psychologically.

## INDICATIONS

### Preoperative preparation

Intestinal absorption problems (tendency to be over-or underweight, linked to the characteristics of the intestinal bacterial flora and the speed of evacuation)  
Haemorrhoids (1<sup>st</sup> and 2<sup>nd</sup> degree)

Scars in the process of healing

Recurring bouts of asthma, bronchitis and pharyngitis

Constipation and alterations of the alvus in general

Skin diseases: acne, mycosis, psoriasis

Weakening of the immune system

Severe headaches

Inflammatory states

Diarrhoea

Intestinal toxemia

Disbiosis

*an age-old cure  
a new technology*



# HYDROGAMMA HCT

*colon hydrotherapy*



## STANDARD ACCESSORIES

Quick couplings for supply

12V power supply unit

Filtering system spare parts

User's manual

10 sets containing speculum and connecting hose,  
in sterile disposable packs

On request: bioresonance kit  
chromo therapy kit  
mineral additive kit  
oxygen and ozone mixing valve

## GENERAL

Ergonomic panel

Compact easily-installed block

Litre counter

Timer

Irrigation pressure monitoring

Bacterial filter

Lighted flow control viewing chamber

## TECHNICAL SPECIFICATIONS

Electrical power supply: 12V DC, from external AC/DC adaptor

Mains water infeed pressure

Automatic pressure reduction

Water outfeed pressure: lower than to blood pressure

Water outfeed temperature: selectable from 32° to 39°C

Classification: CLASS I TYPE B (EN 60601-1)

Marking: CE 0476



# TecnoLaser

Via Venezia 173 • 43122 PARMA - ITALY  
Tel. +39 0521.271102 • Fax +39 0521.272846  
www.tecnolaser.net • e-mail: tecnolsr@tin.it