

TecnoLaser

Manufacturer since 1979

Skin cooling system

Sistema di raffreddamento dell'epidermide



ICE SKIN

ICE SKIN

Skin cooling system

Ideal for laser treatments of dermatologic type, vascular type, laser epilation, for small surgery anaesthesia, multi injections, electro surgery etc.

DOES NOT USE CRYGENIC GASES
ABSENCE OF CONSUMPTION MATERIALS
NO RUNNING COSTS
SIMPLE, EFFICIENT, OF LONG SERVICE LIFE

TECHNICAL SPECIFICATIONS

Minimum temperature:	-10°C
Cooled area	5 cm ²
Colling system	Peltier effect
Power supply	230V - 50/60Hz - 200VA
Dimension	420x250x400 mm
Weight	10 kg
Marking	CE 0476

ICE SKIN

Sistema di raffreddamento dell'epidermide

Ideale per i trattamenti laser di tipo dermatologico, vascolare, epilazione laser e come anestesia per piccoli interventi, multiniezioni, elettrobisturi ecc.

NON UTILIZZA GAS CRIOGENI
ASSENZA DI MATERIALI DI CONSUMO
NESSUN COSTO DI ESERCIZIO
SEMPLICE, EFFICACE, DI LUNGA DURATA

SPECIFICHE TECNICHE

Temperatura minima	-10°C
Area raffreddata	5 cm ²
Sistemi di raffreddamento	effetto Peltier
Alimentazione elettrica	230V - 50/60Hz - 200VA
Dimensioni	420x250x400 mm
Peso	10 kg
Marcatura	CE 0476

TecnoLaser

Via Venezia 173 • 43122 PARMA - ITALY
Tel. +39 0521.271102 • Fax +39 0521.272846
www.tecnolaser.net • e-mail: tecnolaser@tin.it

Handpiece extension for face on request

Prolunga manipolo per il viso a richiesta

