

TecnoLaser

Manufacturer since 1979

From intense pulsed light technology the friendly light of Sirlight for progressive permanent hair removal, photorejuvenation and vascular treatments

Dalla tecnologia della luce super pulsata la luce amica di Sirlight per epilazione progressiva permanente, fotoringiovanimento e trattamento vascolare



SIRIO LIGHT

SIRIO LIGHT From intense pulsed light technology the friend light of Sirlolight for progressive permanent hair removal, photorejuvenation and vascular treatments. Intense pulsed light source equipped with incorporated skin cooling system.

TECHNICAL SPECIFICATIONS

Standard head spectrum	600 - 1200 nm (hair removal)
Optional heads spectrum	on request from 450 to 1200 nm
Fluence	2 - 40 J/cm ²
Pulse duration	5 - 40 ms
Spot dimensions	10 x 50 mm
Power supply	single phase 230V - 50/60Hz
Absorption	6A
Dimensions	370 x 780 x 640 mm
Weight	37 kg
Marking	CE
Optional	new handpiece with interchangeable filters



SIRIO LIGHT Dalla tecnologia della luce super pulsata la luce amica di Sirlolight per epilazione progressiva permanente, fotoringiovanimento e trattamenti vascolari.

Sorgente luce pulsata ad alta intensità con dispositivo di raffreddamento della pelle incorporato.

CARATTERISTICHE TECNICHE

Spettro manipolo standard	600 - 1200 nm (depilazione)
Spettro manipoli optional	su richiesta da 450 a 1200 nm
Fluenza	2 - 40 J/cm ²
Durata impulso	5 - 40 ms
Dimensioni spot	10 x 50 mm
Alimentazione	monofase 230V - 50/60Hz
Assorbimento	6A
Dimensioni	370 x 780 x 640 mm
Peso	37 kg
Marcatura	CE
Optional	nuovo manipolo con filtri intercambiabili

SIRIO LIGHT



 **TecnoLaser**

TECNOLASER srl

Via Venezia 173
43122 PARMA - ITALY
Tel. +39 0521 271102
Fax +39 0521 272846

e-mail: tecnolsr@tin.it

www.tecnolaser.net

www.tecnolaser.de

www.technolaser.eu

www.technolaser.info

www.tecnolasereuropa.info

www.tecnogamma.info

www.applegreen.net

www.applegreen.it

www.applegreen.eu