

TecnoLaser

Manufacturer since 1979

*From intense pulsed light technology
the friendly light of TecnoLight
for progressive permanent
hair removal,
photorejuvenation
and vascular treatments*



*Dalla tecnologia
della luce super pulsata
la luce amica di TecnoLight
per epilazione progressiva
permanente,
fotoringiovanimento
e trattamento vascolare*

TECNOLOGIA

TECNOLIGHT From intense pulsed light technology the friend light of TecnoLight for progressive permanent hair removal, photorejuvenation and vascular treatments.
Intense pulsed light source equipped with incorporated skin cooling system.

TECHNICAL SPECIFICATIONS

Standard head spectrum	600 - 1200 nm (hair removal)
Optional heads spectrum	on request from 450 to 1200 nm
Fluence	2 - 40 J/cm ²
Pulse duration	5 - 40 ms
Spot dimensions	10 x 50 mm
Power supply	single phase 230V - 50/60Hz
Absorption	6A
Dimensions	400 x 400 x 1060 mm
Weight	45 kg
Marking	CE
Optional	new handpiece with interchangeable filters



TECNOLIGHT Dalla tecnologia della luce super pulsata la luce amica di TecnoLight per epilazione progressiva permanente, fotoringiovanimento e trattamenti vascolari.
Sorgente luce pulsata ad alta intensità con dispositivo di raffreddamento della pelle incorporato.

CARATTERISTICHE TECNICHE

Spettro manipolo standard	600 - 1200 nm (depilazione)
Spettro manipoli optional	su richiesta da 450 a 1200 nm
Fluenza	2 - 40 J/cm ²
Durata impulso	5 - 40 ms
Dimensioni spot	10 x 50 mm
Alimentazione	monofase 230V - 50/60Hz
Assorbimento	6A
Dimensioni	400 x 400 x 1060 mm
Peso	45 kg
Marcatura	CE
Optional	nuovo manipolo con filtri intercambiabili

 **TecnoLaser**

Via Venezia 173
43122 PARMA - ITALY
Tel. +39 0521.271102
Fax +39 0521.272846
www.tecnolaser.net
e-mail: tecnolr@tin.it

TECNOLIGHT