

# TecnoLaser

Manufacturer since 1979

## RADIOFREQUENCY AND CRYO The non - invasive new concept in alternative liposculture and lifting

*VELA MAGIE RF looks set to revolutionise the treatment of cellulite and localized fat whilst toning and firming dramatically the treated zone with no risk and side effect. Based on a new multiactive system that combines Radio Frequency Electromagnetic Waves to an exclusive cooling system is a generator of electromagnetic energy at highcomputerized frequency that works towards the membrane of the fat cell to fight , in just few session , adipose cushions and cellulite like anonsurgical alternative to Liposuction and Liposculture.*

*VELA MAGIE RF is furthermore particularly effective in face treatment to contrast skin laxity and ageing , wrinkles and folds.*



# VELA MAGIE RF

It is of fundamental importance to point out that the blood vessel , the lymphatic system and the nerves undergo no alteration in the treatment , thanks to the selectivity of the waves employed which are specific to the different area (body adipose tissue or face).

For the effect of the cold contact with the handpiece the treatment is painless: you feel only a localized heat well within tolerance limits and, at the end, you will feel nothing more than a slight reddening. The result on "orange peel effect" skin , which are visible from the first sessions, are lasting and are Obtained through rehabilitation of the lymphatic vessel s "decompressed" of fat cells which have been reduced in volume.

The lymphatic vessels , freed from compression, very rapidly recover their draining capacity.

**VELAMAGIE RF** is a localized system of silhouette modelling . It may be implemented in a course of a diet to obtain homogeneous slimming with view to overall harmony of the figure , making the skin smooth and well modelled. In the mean time VELAMAGIE RF is a system for alternative lifting to obtain a real anti aging effect , thanks to the stimulation of the collagen production, of the RF emission.

#### CHARACTERISTICS OF THE DEVICE

Power supply	230V 50/60Hz
Absorbed power	240 VA
Dimensions	40 x 40 x 145 (h) cm
Weight	32 kg
Handpiece	n. 1 Monopolar + n.1 Bipolar for face (optional)
Emission power	Adjustable from 5 to 50 W (max)
Source Radiofrequency	with incorporated skin cooling system
Frequency	Adjustable 0,6 -1,5 MHz
Cooling of the handpiece	Adjustable from -5° to +20°C Peltier effect
Cooling of the device	Air-water exchanger.
Emission head diameter	5 cm
Emission head effective	area 20 cm <sup>2</sup>
Classification	I Class BF Type



 **TecnoLaser**

**TECNOLASER srl**  
Via Venezia 173  
43122 PARMA - ITALY  
Tel. +39 0521 271102  
Fax +39 0521 272846  
e-mail: tecnolsr@tin.it

# VELAMAGIE RF